

Foot Surgery



# Operative Technique



Complete plateform
 Low profil plates adapted to each indication
 Presslock<sup>®</sup> Technology





### Table of contents

#### Introduction

| 02 Indications |  |
|----------------|--|
|----------------|--|

02 Contra-indications

### **Design features**

- 03 The Airlock<sup>®</sup> solution
- 04 The Airlock<sup>®</sup> range
- 08 Comprehensive & modular plateform

### **Surgical Technique**

- **16** 1 First MTP Joint Arthrodesis
- 2 Talo-navicular Arthrodesis
- **25** 3 Other indications

### References

- 27 Implants & screws
- 28 Airlock<sup>®</sup> set
- 29 Trial implants
- 31 Distractor set

This document provides technical guidance for the proper usage of the Airlock<sup>®</sup> implant range, however Novastep does not practice medicine and does not recommend this or any other surgical technique. Each surgeon must consider the specific needs of each patient and is responsible for making applicable adjustments and determining and using the appropriate techniques for implanting the device in each given case.

### Introduction

### Indications & Contra-indications

### Indications

The osteosynthesis screw-plate systems are indicated for arthritis (Hallux Rigidus, Osteoarthritis), Hallux valgus and other bone alignment defaults (Hallux Varus, Flatfoot, Cavus foot).

Examples of use:

#### Forefoot:

- . First MTP joint arthrodesis (Hallux-valgus Hallux Rigidus)
- . Revision of first MTP joint arthrodesis
- . Opening or closing Basal osteotomy (Hallux Valgus)
- . Fixation of first metatarsal fractures

#### Midfoot / Rearfoot:

- . Talo-navicular arthrodesis
- . Lapidus arthrodesis
- . Partial or complete Lisfranc arthrodesis
- . Calcaneo-cuboid arthrodesis
- . Evans and Cotton osteotomies
- . Malerba and Dwyer osteotomies
- . Tarsectomy

**Note:** Detailed information on each medical device is provided in the instruction for use. Refer to the instruction for use for a complete list of side effects, warnings, precautions for use, directions for use and contra-indications.



### **Contra-indications**


- . Severe muscular, neurological or vascular deficiency in the extremity concerned.
- . Bone destruction or poor bone quality, likely to impair implant stability.
- . Hypersensitivity to one or more components.

The Airlock<sup>®</sup> plating system is a range of osteosynthesis plates for the forefoot, midfoot and rearfoot, made of TA6V ELI Titanium alloy. It combines a low profile plate design **dedicated to each indication**, a **compression hole**, and a **monoaxial** and **polyaxial** screw system to ensure a stable and rigid fixation.

The Presslock® technology provides a locking compression hole that ensures a stable construct.

### 1 - The Airlock® solution

## Low Profile design > Plates thickness optimized according to the indications to limit subcutaneous discomfort.




### Presslock® innovation - Patent pending

- . Locking compression hole for Ø3.5 mm locking screw.
- . Allows for an additional 1.5 mm compression before locking into a threaded hole, providing strength and stable fixation.



### 2 - The Airlock® range



#### MTP

10°







### Short plates . Thickness: 1.3 mm



#### Standard plates

- · Thickness: 1.3 mm
- Standard compression hole:
   3 mm non-locking screw

#### Long plates


- . Thickness\*: 1.3 mm \*Thickness increases to 1.65 mm
- along the joint line . Standard compression hole:
- 3 mm non-locking screw





1

### LISFRANC



- T plates
- Thickness: 1.5 mm
- Standard compression hole:
   3 mm non-locking screw

### H plates

- Thickness: 1.5 mm
- Standard compression holes:
- 3 mm non-locking screw





### FUSION - Presslock®



### Straight plates

- Thickness: 1.6 mm
- . Presslock<sup>®</sup> compression hole: 3.5 mm locking screw





### H plates

- Thickness: 1.6 mm
- Presslock<sup>®</sup> compression holes:
   3.5 mm locking screws



### LAPIDUS - Presslock®



#### Short plates Long plates . Thickness\*: 1.3 mm

- \*Thickness increases to 1.6 mm in the Presslock® hole area
- Presslock<sup>®</sup> compression hole:
  3.5 mm locking screw



3

4

2

### PLANTAR LAPIDUS - Presslock®



#### Short plates Long plates

#### Long plates

- Thickness: 1.6 mm
  Presslock<sup>®</sup> compression hole:
- 3.5 mm locking screw
- . Transversal screw hole:
- 3.5 mm locking screw



6

6

### UTILITY



#### 6 versions

- 2, 3, 4, 5, 6 & 7 holes
- · Thickness: 1.5 mm
- Standard compression hole:
   3 mm non-locking screw



### MT BASE





MT closing wedge 0 mm Wedge

**MT opening wedge** 3, 4 & 5 mm Wedge Thickness: 1 mm



### MONOAXIAL & POLYAXIAL SYSTEM - Ø 3 & 3.5 mm



\* 2 mm increments.

### COMPRESSIVE SCREWS - Nexis® Ø 4 & 5 mm & PECA®-C Ø 4 mm



\* 2 mm increments from 18 to 50; 5 mm increments from 50 to 60 mm.


\*\* 2 mm increments up to 50 mm, then 5 mm increments.

### 3 - Comprehensive & modular plateform


Used together or separately, the Airlock® instrument trays allow a comprehensive approach to the various indications of foot surgery.







### 3.1 - Airlock® set & optimized instrumentation





#### Windowed drill guides

Direct reading of the required screw length.

#### Drill guides for universal hole:

- . Locking drill guide for drill bits Ø 2 mm & 2.5 mm.
- Polyaxial drill guide for drill bits Ø 2 mm & 2.5 mm.



### Drill guide for


- standard compression hole:
  Compression drill guide for drill bit Ø 2 mm.

### Drill guide for

- Presslock<sup>®</sup> compression locking hole:
- Presslock<sup>®</sup> gold drill guide for drill bits Ø 2 & 2.5 mm.



Drill bits for windowed drill guides



### 3.1.1 - Screw fixation process

The Airlock® Ø 3 and 3.5 mm locking and non-locking screws may be used in all Airlock® plate fixation holes but:

- . Standard compression holes accommodate Ø 3 mm non-locking screws only;
- . Presslock^ ${\ensuremath{\mathbb S}}$  compression locking holes accommodate Ø 3.5 mm locking screws only.

Screw insertion follows an intuitive three-step procedure: drilling, measurement and screw insertion. Each instrument is conveniently organized and color-coded.

#### Color code:

Airlock® screws:

: Instrumentation for Ø 3.5 mm screws

Instrumentation for Ø 3 mm screws

#### Compressive screws:

- : Instrumentation for Nexis® Ø 4 mm screws
- : Instrumentation for Nexis® Ø 5 mm screws
- : Instrumentation for PECA®-C Ø 4 mm screws

Trick: To position a plate, thread 2 locking drill guides in 2 universal holes. Position the plate as desired using the drill guides to manipulate it. Drill the first screw hole with the drill bit of the correct diameter. Leave the drill bit inside to keep the position and drill the second screw hole with a second drill bit. Determine the appropriate screw length by reading the length directly off the windowed drill guide where it matches the calibrated etching on the drill bit or using the depth gauge. Insert the selected screw. Withdraw the drill guide and insert the screw.

Determine the screw length for the other hole before inserting the screw.



In case of uncertainty, screw lengths may be verified by means of the screw length indicator.

### 3.1.1.1 - Universal hole

First, fixate the side of the plate that is opposite to the standard compression hole or Presslock<sup>®</sup> hole, with or without locking. Thread the locking drill guide in one of the threaded holes or position the polyaxial drill guide and drill with the appropriate diameter drill bit. Determine the appropriate screw length by reading the measurement directly off the windowed drill guide or using the depth gauge. Insert the selected screw with the self-retaining screwdriver tip.




#### Instrumentation guidelines



### 3.1.1.2 - Standard compression hole

Start plate fixation opposite the side of the compression hole. The oblong drill guide allows both neutral or compression screw fixation, giving 1.5 mm of additional compression. If no compression is required, use the drill guide in its neutral position.

Note: Compression holes only accomodate the  $\emptyset$  3 mm non-locking screws.




Compression slot generates mechanical compression between the two bone segments.




#### Instrumentation guidelines



Ø 3 mm non-locking screws only



Directly read off the windowed drill guide or use of the depth gauge



Screwdriver tip



#### 3.1.1.3 - Presslock® compression locking hole

Start plate fixation opposite the side of the Presslock® compression hole. The Presslock® drill guide allows neutral or compression screw fixation, giving up to 1.5 mm of compression. To achieve compression locking, position the Presslock® gold drill guide line the Presslock® hole.

If no compression is required, use the Presslock<sup>®</sup> drill guide in its neutral position.

**Note:** Presslock<sup>®</sup> holes, available on the universal Fusion and on the Lapidus and Plantar Lapidus plates, only accommodate Ø 3.5 mm locking screws. Presslock<sup>®</sup> holes are easily identified by their oblong shape and engraved circular arc.



Presslock® slot generates mechanical compression between two bone segments, before subsequent locking in the threaded part of the slot.



Instrumentation guidelines



Ø 3.5 mm locking screw only



Directly read off the windowed drill guide or use of the depth gauge




Screwdriver tip



### 3.1.2 - Concave / Convex reamers

Specific instrumentation, with 3 Concave / Convex reamer sizes (Ø 18 mm, Ø 20 mm, Ø 22 mm) facilitate precise joint surface contouring and positioning. Always be sure to use the concave and convex reamers of the same diameter.



### 3.1.3 - Plate bender procedures

Most of the time, bending is not necessary. In some rare cases plate benders may be required. The following guidelines must be considered:

- . Bend the plate only in one direction.
- . Never reverse-bend a plate.
- . Always ensure that the threaded holes of a plate are not compromised during bending.
- . It is not recommended to bend at the plate extremities.







Note: The following plates must not be bent in order to avoid damaging the central hole: . MTP Long . Plantar Lapidus

#### 3.2 - Distractors & compressor set

#### **5** Distractors & compressor



### 3.2.1 - Opening Wedge distractor




- (1) Anatomical wedge-shaped opening Made possible by the ball-and-socket joint on the dedicated threaded K-Wires
- (2) Opening Wedge threaded K-wires Allow a stable opening and a distraction close to the bone
- (3) Easy measurement of the opening After K-Wires insertion adjust the initial position to 0mm The sliding graduated ruler allows 0-12 mm sizing increments
- (4) Simple and precise adjustment Of the opening, fixed by a screwed wheel

### Example of the use of the Opening Wedge distractor on an Evans osteotomy

Placement of the wires and distractor at the bone cut

Setting the initial position to 0 mm on the ruler before distraction: Turn the dial (1) to allow the movement of the ruler (2).



Performing the distraction with a wedge-shaped anatomical opening



Placement of a bone graft and an Airlock® Fusion plate



### 1 - First MTP joint arthrodesis

#### MTP plating system benefits

3 sizes available in Left & Right option: Short, Medium & Long

41 mm





- Threaded holes for Ø 3 mm & Ø 3.5 mm locking or nonlocking screws
- The anatomic plates are designed with 0° dorsiflexion, delivering 15° of metatarsophalangeal dorsiflexion while preserving 10° anatomical phalangeal valgus

Low-Profile plate design reduces soft tissue irritation around the MTP joint


Valgus 10°

53 mm

#### Thickness: 1.3 mm

34 mm

MTP Long: for improved strength, the plate increases in thickness from 1.3 mm to 1.6 mm along the joint line



### 1 - Incision & exposure

A medial incision is most commonly used for First MTP joint exposure. A dorsal approach could be also considered. It is recommended to identify and to protect the dorsal collateral nerve to avoid any risk of damage during opening or closing steps. An exostectomy is performed with an oscillating saw and a large circumferential arthrolysis is performed to expose the entire joint area. Osteophytes are completely resected.

### 2 - Metatarsal & phalangeal preparation

Care should be taken to protect skin and soft tissue during the joint surface preparation. Two options may be considered: Flat Cut Technique or Cup & Cone Technique.

- The Cup & Cone Technique allows an easier and more precise adjustment that preserves bone stock, but requires more exposure. Adapt the cut with consideration for first ray length and overall bone quality (cancellous and sclerotic bones).
- . In case of shortening, flat cuts are recommended.
- . In case of poor quality bones or osteoporotic bones, gouge forceps preparation is preferred.
- In case of sclerotic bone, the Cup & Cone Technique is recommended. To facilitate fusion, weaken the counteropposing surfaces with gouge forceps, oscillating saw or bone scraper prior to application of the Cup & Cone reamers.


Displace the phalanx plantarly, exposing the metatarsal head (it is recommended to start the metatarsal preparation first to enable proper exposure of the phalanx).

Using a power drill, place a Ø 1.6 mm K-Wire through the center of the metatarsal head and into the diaphysis of the metatarsal. Utilize the largest reamer size to start the metatarsal reaming process.

Reaming of the phalanx is performed in a similar fashion to the metatarsal head.

**Note:** Joint surface may be prepared by performing perforations with a drill bit or K-wire.







### 3 - Temporary fixation

Provisionally stabilize the joint by inserting a 1.6 mm K-Wire from the dorsal medial aspect of the first metatarsal to the dorsal lateral cortex of the first phalanx.

Check the correct position (approximately 15° of dorsiflexion) using the support plate, located inside the lid of the instrument tray as pictured below:

- . Dorsal Aspect: Hallux is parallel to the second toe and nail is parallel to the ground.
- Medial / Lateral Aspect: With the patient's heel resting on the plate, the pad of the big toe should be slightly elevated (<5 mm). The great toe needs to have the ability to stay in contact with the floor.



**Medial / Lateral Aspect:** Difference between the support plate and the pulp of the big toe < 5 mm.




### 4 - Plantar screw insertion

**4.1** - Place the trial implant using the locking drill guide screwed in the proximal hole. Plantar screw guide marks on the surface of the trial implant depict the recommended orientation for insertion of the Ø 1.4 K-Wire (approximately  $35^{\circ}$  relative to the metatarsal axis). Identify the appropriate plate reference according to the trial implants.



**4.3 - Option 1: Self-drilling screw:** Use the T10 screwdriver tip to insert the Ø 4 mm compressive screw manually (using a driver handle) or with a power tool. Check to ensure proper stability at the osteotomy site.

4.2 - Determine screw length using the measuring gauge.



**4.4 - Option 2: Pre-drilling and Countersink:** Prepare bone housing using the dedicated  $\emptyset$  2.7 mm Nexis<sup>®</sup> drill bit and  $\emptyset$  3.7 countersinking reamer. Insert the screw with the T10 screwdriver tip.



**Optional:** Replace the Nexis<sup>®</sup> 4 mm screw with the PECA<sup>®</sup> Compressive 4 mm screw to allow maximization of cortical anchorage and preservation of soft tissues thanks to the bevelled head. In this case, be sure to use the associated instrumentation: the Exact-T<sup>®</sup>10 screwdriver tip and, if necessary, the PECA<sup>®</sup> 4 Ø 3.2 mm AO drill bit and the Nexis<sup>®</sup> / PECA<sup>®</sup> Ø 3.7 mm countersink.






### 5 - Plate positioning

**5.1 - Plate housing preparation:** If necessary, flatten the dorsal surface using a oscillating saw or gouge forceps.


**5.2 - Plate positioning:** With the joint now stabilized, the plate should be placed over the joint and positioned according to patient's anatomy. When the proper orientation is determined, insert the spheric positioning pins to secure the plate over the bone as illustrated below.



**5.3 - Distal screw insertion:** The steps for inserting Airlock<sup>®</sup> screws and the use of associated instruments are specified in the Introduction (paragraphs 3.1.1.1 and 3.1.1.2).

Prepare the most distal screw hole using the drill bit  $\emptyset$  2.0 mm and locking drill guide (for locking screws) or the polyaxial drill guide (for non-locking screws).

Determine the appropriate screw length by directly reading off the windowed drill guide or using the depth gauge. Insert the selected screw.




**Note:** It is recommended to first insert a non-locking screw prior to introducing locking screws to position the plate flush with respect to the cortical surface. It is also recommended to perform distal fixation prior to inserting the proximal screws and always prior to using the proximal compression hole.


### 6 - Proximal screw insertion

Using the compressive drill guide, drill the compression screw hole for non-locking screw.

Each proximal wire must be removed prior to compressing the joint.

If compression is not required, use the compressive drill guide in its neutral position (see Paragraph 3.1.1.2 of the Introduction). Determine the appropriate screw length and insert the screw until full compression is achieved.





### 7 - Screw insertion

Insert remaining screws and check the stability of the assembly.



### 2 - Talo-navicular arthrodesis

### 1 - Incision & exposure

A dorsal incision is most commonly used for TN joint exposure. A medial incision can be made following the surgeon's preference.

Once the joint is exposed, position a closed arm distractor on the talus and navicular bones. Sterile threaded wires are available for use with the distractor.

Distract the joint and remove articular cartilage using curette, rongeur or small osteotome.

A K-wire can be placed across the talo-navicular joint to stabilize the joint.



### 2 - Trial implants

Use the trial Presslock<sup>®</sup> universal Fusion plates to determine the appropriate shape and size.


Depending on the joint anatomy, 1 or 2 Straight plates or 1 H plate can be used, all available in Short, Medium or Long sizes.



### 3 - Inter-fragmentary compression screw insertion

Insert a Ø 1.4 mm K-wire from the navicular to the talus and determine the screw length using the depth gauge.





**Option 1: Self-drilling screw:** Use the T10 screwdriver tip to insert the  $\emptyset$  4 mm compressive screw manually (using a driver handle) or with a power tool. Check to ensure proper stability at the osteotomy site.

**Option 2: Pre-drilling and Countersink:** Prepare bone housing using the dedicated  $\emptyset$  2.7 mm Nexis<sup>®</sup> drill bit and  $\emptyset$  3.7 countersinking reamer. Insert the screw with the T10 screwdriver tip.



**Optional:** Replace the Nexis<sup>®</sup> 4 mm screw with the PECA<sup>®</sup> Compressive 4 mm screw to allow maximization of cortical anchorage and preservation of soft tissues thanks to the beveled head. In this case, be sure to use the associated instrumentation: the Exact-T<sup>®</sup>10 screwdriver tip and, if necessary, the PECA<sup>®</sup> 4 Ø 3.2 mm AO drill-bit and the Nexis<sup>®</sup> / PECA<sup>®</sup>-C countersink Ø 3.7 mm.

#### 4 - Plate positioning

Position the Presslock<sup>®</sup> universal Fusion plate according to the patient's anatomy. Secure the plate with temporary fixation pin.

**Screw insertion:** The steps for inserting Airlock<sup>®</sup> screws and the use of associated instruments are specified in the Introduction - paragraphs 3.1.1.1 and 3.1.1.3.

**4.1 - Universal hole:** Ø **3.5 mm screw insertion:** Prepare the screw hole using either the locking drill guide for locking screws or the polyaxial drill guide for non-locking screws.

Determine the appropriate screw length by either directly reading through the windowed drill guide or using the depth gauge after removing the drill guide, and insert the appropriate screw.



#### **4.2 - Presslock® hole:** Ø **3.5 mm locking screw insertion:** Position the Presslock® drill guide into the compression locking Presslock® hole. Prepare the screw hole.

Determine the appropriate screw length by either directly reading through the windowed drill guide or by using the depth gauge after removing the drill guide and insert the locking screw into the locking compression slot.

Repeat these steps for the second  $\mathsf{Presslock}^{\circledast}$  hole if a H-plate has been used.

Check the construct stability and confirm placement using fluoroscopy.





4.3 - Final positioning:



H plate final positioning



Straight plates final positioning

### 3 - Other indications

Regardless of the Airlock® plate used for different indications, follow the same steps for inserting Airlock® screws and the use of associated instruments specified in the Introduction.

### Examples of Airlock<sup>®</sup> plating system:

### 1 - Lisfranc arthrodesis



Lisfranc H-plate



Lisfranc T-plate

### 2 - Lapidus arthrodesis



Presslock® Lapidus plate




Presslock® plantar Lapidus plate



Presslock<sup>®</sup> Fusion straight plates

3 - Closing wedge osteotomy



MT closing wedge plate - Wedge 0 mm

### 4 - Calcaneocuboïd arthrodesis



Presslock<sup>®</sup> Fusion H plate

### 5 - Naviculocuneiform arthrodesis



Presslock<sup>®</sup> Fusion H plate

### 6 - Evans osteotomy



Presslock<sup>®</sup> Fusion H plate

### 7 - Malerba osteotomy



Presslock® Fusion straight plate

### 8 - Dwyer osteotomy



Presslock® Fusion straight plate

### 1 - Implants & screws

#### MTP

| Reference      | Designation           |
|----------------|-----------------------|
| PL010134 / 234 | Short - Right / Left  |
| PL010140 / 240 | Medium - Right / Left |
| PL010152 / 252 | Long - Right / Left   |

#### Fusion - Presslock®

| Reference            | Designation                      |
|----------------------|----------------------------------|
| PL040117 / 120 / 123 | Straigth – Short / Medium / Long |
| PL040217 / 220 / 223 | H – Short / Medium / Long        |

### Lapidus - Presslock®

| Reference     | Designation          |
|---------------|----------------------|
| PL030301 / 02 | Short – Right / Left |
| PL030401 / 02 | Long – Right / Left  |

### Plantar Lapidus - Presslock®

| Reference      | Designation          |
|----------------|----------------------|
| PL080101 / 102 | Short – Right / Left |
| PL080201 / 202 | Long – Right / Left  |

#### Airlock<sup>®</sup> screws Ø 3 & Ø 3.5

| Locking  |          | Non l    | ocking   |
|----------|----------|----------|----------|
| Ø 3.0 mm | SP0130XX | Ø 3.0 mm | SP0230XX |
| Ø 3.5 mm | SP0135YY | Ø 3.5 mm | SP0235YY |

With XX from 10 to 30 in 2 mm increments. With YY from 10 to 40 in 2 mm increments.

#### **Compressive screws**

| Screws        | Length       | Reference |
|---------------|--------------|-----------|
| Nexis® Ø 4.0  | 18 to 60 mm  | SC0500XX  |
| Nexis® Ø 5.0  | 30 to 100 mm | SC060YYY  |
| PECA®-C Ø 4.0 | 18 to 60 mm  | PS0501XX  |

With XX from 18 to 60 in 2 mm increments up to 50 and 5 mm increments up to 60.

With YYY from 030 to 050 in 2 mm increments up to 50 and 5 mm up to 100.

### Screws color code

| Color | Screws                        |
|-------|-------------------------------|
| •     | Airlock <sup>®</sup> Ø 3.0 mm |
| •     | Airlock® Ø 3.5 mm             |
| •     | Nexis® Ø 4.0 mm               |
| 0     | Nexis® Ø 5.0 mm               |
| •     | PECA®-C Ø 4.0 mm              |

### Lisfranc

| Reference | Designation |
|-----------|-------------|
| PL050101  | T - Short   |
| PL050102  | T - Long    |
| PL050201  | H - Small   |
| PL050202  | H - Medium  |
| PL050203  | H - Large   |

#### Utility

| Reference | Designation        |
|-----------|--------------------|
| PL040016  | 2 Holes, Length 16 |
| PL040022  | 3 Holes, Length 22 |
| PL040028  | 4 Holes, Length 28 |
| PL040034  | 5 Holes, Length 34 |
| PL040040  | 6 Holes, Length 40 |
| PL040046  | 7 Holes, Length 46 |

### MT Base


| Reference      | Designation                     |
|----------------|---------------------------------|
| PL020100 / 200 | MT closing Wedge - Right / Left |
| PL020103 / 203 | MT open Wedge 3 - Right / Left  |
| PL020104 / 204 | MT open Wedge 4 - Right / Left  |
| PL020105 / 205 | MT open Wedge 5 - Right / Left  |

### K-Wires & positioning pins

| Reference       | Airlock <sup>®</sup> tray     |
|-----------------|-------------------------------|
| XPP01003        | Spherical positioning pin     |
| 33-T10-R-14-150 | K-Wire Ø 1.4 Lg 100           |
| 33-T10-R-16-180 | K-Wire Ø 1.6 Lg 180           |
| Reference       | Distractors & Compressor trav |

| Herenenee | Biolidotoro di Compresser il dy                            |
|-----------|--|
| SKW04001  | Threaded K-Wire Ø 2.5 - opening wedge distractor - sterile |
| SKW05003  | Threaded K-Wire TR-RD Ø 1,6 lg 150 - sterile               |
| SKW05004  | Threaded K-Wire TR-RD Ø 2,5 lg 200 - sterile               |

### 2 - Airlock<sup>®</sup> set



### Empty tray / lid

| Ref          | Description        |   |
|--------------|--------------------|---|
| ACC1016P0001 | Tray               | Α |
| ACC1016P0003 | Lid                | В |
| ACC1002P0008 | Foot support plate | С |

### **Universal instruments**

| Ref          | Description        | Image |
|--------------|--------------------|-------|
| XHA01001     | AO handle          |       |
| XGA01002     | Depth gauge        |       |
| ACC1016P0004 | K-Wires holder*    |       |
| XKW01002     | Cleaning pin Ø 1.4 | 0     |

\* K-Wire (33-T10-R-14-150): Ø 1.4 for Nexis® Ø 4.0 mm. (33-T10-R-16-180): Ø 1.6 mm for Airlock® reamers.

### Airlock® instrumentation

| Ref      | Description                                   | Image     |
|----------|---|-----------|
| XSD02002 | T8 AO screwdriver tip                         |           |
| XPP01003 | Spheric positioning pin                       |           |
| XGA01003 | Screw length indicator                        |           |
| XDG01020 | Polyaxial drill guide for screw Ø 3 and Ø 3.5 | - <u></u> |
| XDG01021 | Locking drill guide for screw Ø 3 and Ø 3.5   |           |
| XDG01022 | Oblong drill guide for screw Ø 3              |           |
| XDG01023 | Presslock® - Gold drill guide for screw Ø 3.5 |           |
| XDB01021 | Drill bit Ø 2                                 |           |
| XDB01022 | Drill bit Ø 2.5                               |           |
| XMS01004 | Drill guide holder                            | 5 mmm     |

### **Reamers & plate benders**

| Ref      | Description         | Image  |
|----------|---------------------|--|
| XMS01010 | Plate benders       |  |
| XRE01010 | Convex reamer Ø 18  | · · · · · · · · · · · · · · · · · · ·  |
| XRE01005 | Convex reamer Ø 20  | ·  |
| XRE01006 | Convex reamer Ø 22  |  |
| XRE01011 | Concave reamer Ø 18 | The second secon |
| XRE01003 | Concave reamer Ø 20 |  |
| XRE01004 | Concave reamer Ø 22 |  |

### Nexis<sup>®</sup> & PECA<sup>®</sup>-C instrumentation

| Ref      | Description   | Image  |
|----------|---|--|
| XSD03001 | Snap-off AO screwdriver tip                           | XSD00001 🛷   |
| XSD04001 | T10 AO screwdriver tip                                | x306401 🛷 📑  |
| XRE01007 | Nexis® / PECA®-C countersink Ø 3.7                    |  |
| XDB01007 | Cannulated drill bit Ø 2.7                            |  |
| XDG01009 | Double drill guide for screw Ø 4                      | and the first state of the stat |
| XGA01009 | Nexis <sup>®</sup> / PECA <sup>®</sup> - Ruler Lg 150 | 0 0 0 0 4 stantes<br>6 0 6 6 g f terms   |
| XSD05001 | T20 AO screwdriver tip                                | Teleo X10c10c1 🛩   |
| XRE01008 | Countersink Ø 4.9                                     | XRESIDS 🛷  |
| XDB01009 | Cannulated drill bit Ø 3.2                            |  |
| XDG01015 | Double drill guide for screw Ø 4                      | And the free free to be  |
| XGA01007 | Ruler Lg 180/200                                      |  |
| CKW02005 | K-wire Ø 1.4 Lg 150 TR/RD CrCo                        |  |
| XSD04004 | Exact-T®10 AO screwdriver tip                         | Exect 110  |
| XDB01023 | Drill bit Ø 3.2                                       |  |

### 2.1 - Trial implants

### **MTP** trial plates

| Ref          | Description                  | Image           |
|--------------|------------------------------|-----------------|
| ACC1006P0009 | MTP trial plate holder       |                 |
| XTI01301     | MTP short trial plate left   | COM             |
| XTI01302     | MTP short trial plate right  |                 |
| XTI01001     | MTP medium trial plate left  | Commente        |
| XTI01002     | MTP medium trial plate right | OR              |
| XTI01401     | MTP long trial plate left    | dig Dig no      |
| XTI01402     | MTP long trial plate right   | <b>Citagino</b> |

### **Universal Fusion trial plates**

| Ref          | Description                        | Image      |
|--------------|------------------------------------|------------|
| ACC1006P0007 | Fusion trial plate holder          |            |
| XTI04117     | Fusion straight short trial plate  | (C))(C)    |
| XTI04120     | Fusion straight medium trial plate | Canado     |
| XTI04123     | Fusion straight long trial plate   |            |
| XTI04217     | Fusion H short trial plate         | 2          |
| XTI04220     | Fusion H medium trial plate        | <b>H</b>   |
| XTI04223     | Fusion H long trial plate          | <b>2-6</b> |

### Lisfranc trial plates

| Ref          | Description                   | Image                   |
|--------------|-------------------------------|-------------------------|
| ACC1006P0012 | Lisfranc trial plate holder   |                         |
| XTI05010     | Lisfranc H small trial plate  | <b></b>                 |
| XTI05020     | Lisfranc H medium trial plate | <b></b>                 |
| XTI05030     | Lisfranc H large trial plate  | 11                      |
| XTI05040     | Lisfranc T short trial plate  | 600000G                 |
| XTI05050     | Lisfranc T long trial plate   | accounting and a second |

### Lapidus & plantar lapidus trial plates

| Ref          | Description                             | Image                                   |
|--------------|---|---|
| ACC1006P0011 | Lapidus trial plate holder              |   |
| XTI03301     | Lapidus short trial plate right         | • · · · · · · · · · · · · · · · · · · · |
| XTI03302     | Lapidus short trial plate left          |   |
| XTI03401     | Lapidus long trial plate right          | C A A A A A A A A A A A A A A A A A A A |
| XTI03402     | Lapidus long trial plate left           | Q BEERTO,                               |
| XTI08101     | Plantar Lapidus short trial plate right |   |
| XTI08102     | Plantar Lapidus short trial plate left  | ······································  |
| XTI08201     | Plantar Lapidus long trial plate right  | •10 k line *                            |
| XTI08202     | Plantar Lapidus long trial plate left   |   |


### Utility trial plates

| Ref          | Description                 | Image     |
|--------------|-----------------------------|-----------|
| ACC1006P0010 | Utility trial plate holder  |           |
| XTI04016     | Utility 2 holes trial plate | CINED     |
| XTI04022     | Utility 3 holes trial plate |           |
| XTI04028     | Utility 4 holes trial plate | COMMEND   |
| XTI04034     | Utility 5 holes trial plate |           |
| XTI04040     | Utility 6 holes trial plate | GCG31200  |
| XTI04046     | Utility 7 holes trial plate | ODC0==200 |

### MT Base trial plates - optional

| Ref          | Description                        | Image  |
|--------------|------------------------------------|--|
| ACC1006P0008 | Basal MT trial plate holder        |  |
| XTI02010     | MT Closing wedge trial plate left  | diamon and the second s |
| XTI02013     | MT Open wedge 3 trial plate left   | Amo  |
| XTI02014     | MT Open wedge 4 trial plate left   | dence o  |
| XTI02015     | MT Open wedge 5 trial plate left   | diamon and the second  |
| XTI02020     | MT Closing wedge trial plate right | cond,  |
| XTI02023     | MT Open wedge 3 trial plate right  | cond,  |
| XTI02024     | MT Open wedge 4 trial plate right  | crange   |
| XTI02025     | MT Open wedge 5 trial plate right  | could  |

### 3 - Distractor set



| Ref          | Description  | Image           |
|--------------|--|-----------------|
| ACC1016P0002 | Distractor tray  |                 |
| ACC1016P0005 | Distractor K-Wires holder                                  |                 |
| CKW04001     | Threaded K-Wire Ø 2.5 - opening wedge distractor - sterile |                 |
| SKW05003     | Threaded K-Wire Ø 1.4 lg 150 TR / RD - sterile             |                 |
| SKW05004     | Threaded K-Wire Ø 1.6 lg 150 TR / RD - sterile             |                 |
| XFP01006     | Closed arms distractor                                     | $\triangleleft$ |
| XFP01008     | Open arms distractor                                       | $\triangleleft$ |
| XFP01011     | Opening wedge distractor                                   | $\triangleleft$ |
| XFP01012     | Compressive forceps  | 2               |

### **Notes**

### **Notes**

### **Notes**



#### Please Note:

Carefully read the enclosed Instructions For Use (IFU) and all packaging label information. Devices: Implants: Class IIb-CE1639 / Instruments: Class I / Class IIa-CE1639.

Novastep:
 2, Allée Jacques Frimot - 35000 RENNES - France
 Tel: + 33 (0) 2 99 33 86 50 / Fax: + 33 (0) 9 70 29 18 95
 contact@novastep-ortho.com / www.int.novastep.life

Reference: Air-ST-Ed6-02-24-EN