

arcad[®]




BROCHURE

FOREFOOT
MIDFOOT
REARFOOT

Osteosynthesis compressive staples

8 / 10 / 15 / 18 / 20 / 25 mm



- . No temperature management
- . Fast, intuitive insertion
- . Dynamic compression
- . Mono-cortical or bi-cortical anchorage


*Creating
Better
Together™*



Made of Nitinol, osteosynthesis compressive staples are indicated for the correction of bone deformities (Hallux valgus, Lesser ray deformities, Hammertoe deformity, Coalitions, Rheumatoid foot deformity, Cavovarus foot, Flatfoot), and treatment of foot joint stiffness (Hallux rigidus, Foot arthritis).

Designed to facilitate fast and easy fixation, memory metal staples provide dynamic compression at the osteotomy and arthrodesis site.

Arcad® Solution !



Features

Shape Memory Alloy for immediate and dynamic compression.

No temperature management.

Compressive strength is generated by mechanical stress by means of Arcad® insertion forceps.

Low profile / atraumatic contours.

Superelastic Effect

Nitinol is an alloy made of 50% Nickel and 50% Titanium. Nitinol implants have the ability to return to their original shape immediately after unloading, inducing immediate and dynamic inter-fragmentary compression.



Range


Size range allows for precise fitting and anatomic compatibility.

- Width options : 8 / 10 / 15 / 18 / 20 / 25 mm.

- Mono-Cortical or Bi-Cortical Anchorage.

Raw Material: Nitinol.

Sizing decision-making support - Sizes example by procedure*



Forefoot :

1. Akin Osteotomy	08/10
2. Bunion	10
3. MTP Arthrodesis	15/18

Midfoot :

4. Naviculocuneiform Arthrodesis	15/18
5. Lisfranc Arthrodesis	15/18
6. Lapidus Arthrodesis	15/18

Rearfoot :

7. Talo-Navicular Arthrodesis	20/25
8. Calcaneo-Cuboidal Arthrodesis	20/25

MTP Arthrodesis



Akin Osteotomy



Lisfranc Arthrodesis




*Size choice example to be adjusted at the surgeon's discretion according to the patient anatomy.

Minimal Instrumentation, Simplified Technique


Intuitive Instrumentation

Fast insertion : Drill, Measure, Deploy and Insert. Repositioning and per-operative removal facilitated.

Arcad® Module



Arcad® Module
15/18/20/25




RearfootCOMPLETE

Modular Solution


The Arcad® instrumentation is delivered in specific, intuitively configured modules for forefoot, midfoot and rearfoot applications by means of the EXACT and COMPLETE sterilization tray platforms, which include the Lync intramedullary implants, Arcad® staples, Nexis screws and Percutaneous Reamers.

ForefootEXACT



Arcad® Module
8/10/15/18

Arcad® Module
8/10/15/18




ForefootCOMPLETE



Arcad® Forceps : One Click Expansion !

Arcad® Forceps deliver a specific compression strength for each implant size. The implantation forceps for the Arcad® Staples feature Novastep's unique Lock and Load design that permits convenient, reproducible deployment and rapid insertion.



References



Arcad[®] - Compressive staples

References	Description	Bridge Width (mm)	Leg Length (mm)	Wire Dimensions (mm)
CS020808	Compressive Staple 08-08-08	8	8-8	1,2 x 1,5
CS010909	Compressive Staple 10-09-09	10	9-9	1,2 x 1,5
CS011214	Compressive Staple 10-12-14	10	12-14	1,2 x 1,5
CS011416	Compressive Staple 10-14-16	10	14-16	1,2 x 1,5
CS011518	Compressive Staple 10-15-18	10	15-18	1,2 x 1,5
CS031212	Compressive Staple 15-12-12	15	12-12	1,3 x 1,8
CS031414	Compressive Staple 15-14-14	15	14-14	1,3 x 1,8
CS041212	Compressive Staple 18-12-12	18	12-12	1,3 x 1,8
CS041414	Compressive Staple 18-14-14	18	14-14	1,3 x 1,8
CS041416	Compressive Staple 18-14-16	18	14-16	1,3 x 1,8
CS041618	Compressive Staple 18-16-18	18	16-18	1,3 x 1,8
CS041818	Compressive Staple 18-18-18	18	18-18	1,3 x 1,8
CS051818	Compressive Staple 20-18-18	20	18-18	2,5 x 1,6
CS062020	Compressive Staple 25-20-20	25	20-20	2,5 x 1,6

In addition to the Arcad[®] range

A straight and angulated stainless steel staple range for variation osteotomies



Static staples

References	Description
FS010008	Static Staple 26° - 8 mm
FS010010	Static Staple 26° - 10 mm
FS010108	Static Staple 90° - 8 mm
FS010110	Static Staple 90° - 10 mm

Please note:

Carefully read the enclosed Instructions For Use (IFU) and all packaging label information. Devices: Implants: Class IIb-CE1639 / Instruments: Class I / Class Ir-CE1639 / Class IIa-CE1639.

Novastep:

2, Allée Jacques Frimot - 35000 RENNES - France
 Tel: + 33 (0) 2 99 33 86 50 / Fax: + 33 (0) 9 70 29 18 95
 contact@novastep-ortho.com / www.int.novastep.life

Reference: Arc-L-Ed6-07-24-EN