

NEXIS[®] FOREFOOT


SCREW SYSTEMS



KEY FEATURES

Wide range of titanium alloy compressive screws, designed with self-drilling and self-tapping properties to provide a quick and effective solution intended for the fixation of osteotomy and arthrodesis of the forefoot.

FOREFOOT SOLUTION



INDICATIONS

The osteosynthesis screws are indicated for arthrosis, hallux valgus and other bone alignment defaults (pes cavus, flatfoot, malalignment secondary to previous trauma).



EXAMPLE OF USE

- Ø2mm: Weil osteotomy
- Ø2.3mm and Ø2.9mm: P1 osteotomy, Scarf & Chevron
- Ø4mm: MTP & Cuneo-Metatarsal arthrodesis

MATERIALS


Titanium Alloy TA6V ELI anodized

MODULAR INSTRUMENTATION

A modular concept allowing customization of the ancillary according to the preferences of the practitioner or the indication.

There are two versions available:

- **ForefootCOMPLETE:** possibility of accommodating Nexis® screw, Arcad® staple and Lync® intramedullary implant modules in a complete tray.
- **ForefootEXACT:** possibility of interchanging the modules Nexis®, Arcad® staple and Lync® intramedullary implant for space saving in a half-tray.



FOREFOOTCOMPLETE




FOREFOOTEXACT

1. Nexis® screw module
2. Lync® intramedullary implant module
3. Arcad® staple module

DEDICATED INSTRUMENTATION

DEDICATED INSTRUMENTATION

Dedicated instrumentation for each osteotomy: Weil, Scarf, Chevron, M1.



SNAP-OFF AO SCREWDRIVER TIP


Snap-off screws compatible for use with wire driver or for manual fixation with specific snap-off screwdriver tip.



**SPOON FOR
WEIL OSTEOTOMY**
(XMS01007)*




**SCARF
FORCEPS**
(XFP01001)*



**CHEVRON
FORCEPS**
(XFP01004)*

COLOR CODED


Color coded equipment permits at-a-glance instrument identification.




● Ø2MM



○ Ø2.3MM



● Ø2.9MM



● Ø4MM

ORDERING INFORMATION

SNAP-OFF SCREWS

LENGTH (mm)	Ø2 mm
11	SC040011
12	SC040012
13	SC040013
14	SC040014
15	SC040015
16	SC040016

K-WIRES

PART NO.	DESCRIPTION
-	K-WIRE Ø0,9 LG 80 TR/RD ⁽¹⁾
-	K-WIRE Ø 1,0 LG 80 TR/RD ⁽²⁾
-	K-WIRE Ø 1,4 LG 100 TR/RD ⁽³⁾

⁽¹⁾K-wire supplied separately - Medetechnik® K-wire (33-T10-R-09-080) or Novastep® K-wire (CKW01009) are available depending on your market.

⁽²⁾K-wire supplied separately - Medetechnik® K-wire (33-T10-R-10-080) or Novastep® K-wire (CKW01016) are available depending on your market.

⁽³⁾K-wire supplied separately - Medetechnik® K-wire (33-T10-R-14-100) or Novastep® K-wire (CKW01002) are available depending on your market.

COMPRESSIVE SCREWS

LENGTH (mm)	Ø2.3 mm	Ø2.9 mm	Ø4 mm
10	SC010010	SC020010	-
12	SC010012	SC020012	-
14	SC010014	SC020014	-
16	SC010016	SC020016	-
18	SC010018	SC020018	SC050018
20	SC010020	SC020020	SC050020
22	SC010022	SC020022	SC050022
24	SC010024	SC020024	SC050024
26	SC010026	SC020026	SC050026
28	SC010028	SC020028	SC050028
30	SC010030	SC020030	SC050030
32	-	SC020032	SC050032
34	-	SC020034	SC050034
36	-	-	SC050036
38	-	-	SC050038
40	-	-	SC050040
42	-	-	SC050042
44	-	-	SC050044
46	-	-	SC050046
48	-	-	SC050048
50	-	-	SC050050
55	-	-	SC050055
60	-	-	SC050060

NOTES

enovis.

T +33 (0) 2 99 33 86 50 F + 33 (0) 9 70 29 18 95

Legal manufacturer: Novastep® S.A.S
2 Allée Jacques Frimot | 35000 Rennes | France
contact-nsp@enovis.com
www.int.novastep.life

Copyright © 2024 Enovis Foot and Ankle

Carefully read the enclosed Instructions For Use (IFU) and all packaging label information. Devices: Implants: Class IIb-CE1639 / Instruments: Class I / Class Ir-CE1639 / Class IIa-CE1639.

REFERENCE: NEX-FORE-L-ED7-06-25-EN