

**nexis** Nexis Module 2.0/2.3/2.9/3.5

Part#	Description	QTY
ACC1001P0006	Module Nexis	1
<b>1</b> XSD01001	Nexis AO screwdriver shaft T7	1
<b>2</b> XSD02001	Nexis AO screwdriver shaft T8	1
<b>3</b> XSD03001	Nexis snap-off screwdriver	1
<b>4</b> XRE01001	Nexis countersink reamer	1
<b>5</b> XDB02001	Solid drill bit Ø1.75	1
<b>6</b> XDB01001	Canulated drill bit Ø1.75	1


**arcad** Arcad Forefoot Module

Part#	Description	QTY
ACC1001P0003	Arcad caddy	1
<b>1</b> XFP03001	Arcad 10 forceps	1
<b>2</b> XDG01001	Arcad 10 drill guide	1
<b>3</b> XDB01008	Arcad drill bit Ø2	1
<b>4</b> XFP02001	Impactor staples angled	1
<b>5</b> XMS01002	Staple impactor	1
<b>6</b> XPP01001	Positioning pin Ø2	2
<b>7</b> XFP05001	Impactor staple straight	1

**Universal Instruments**

Part#	Description	QTY
ACC1001P0001	Forefoot tray	
ACC1001P0002	Forefoot Lid	
<b>1</b> XFP01001	Scarf Forceps**	1
<b>2</b> XFP01004	Chevron Forceps**	1
<b>3</b> XHA01001	AO Handle	1
<b>4</b> XKW01001	Cleaning wire	1
<b>5</b> XGA01001	Graduated ruler	1
<b>6</b> XGA01002	Depth gauge	1
<b>7</b> XMS01001	K-wire tube*	1

\*Ref K-wire: 33-T10-R-09-080: K-wire Ø 0.9 mm lg 80 mmTr/Rd • 33-T10-R-10-080: K-wire Ø 1.0 mm lg 80 mmTr/Rd • 33-T10-R-14-100: K-wire Ø 1.4 mm lg 100 mmTr/Rd  
 \*\*Optional




Optional: Nexis 4 Module (page 2)

**lync** Lync Module

Part#	Description	QTY
ACC1001P0004	Lync caddy	1
<b>1</b> XFP04001	Lync forceps	1
<b>2</b> XDB01003	Drill bit Ø2.3	1
<b>3</b> XRE01002	Lync surfacing reamer	1
<b>4</b> XRA01002	Lync distal rasp	1

**nexis** Nexis 4.0 Module



**nexis** Nexis 4.0 Module  
**Part#**

	<b>Part#</b>	<b>Description</b>	<b>QTY</b>
	ACC1002P0004	Nexis Ø4 module	1
<b>1</b>	XSD04001	Nexis AO T10 screwdriver tip for 4.0mm screws	1
<b>2</b>	XDG01009	Nexis double drill guide 4.0mm screws	1
<b>3</b>	XDB01007	Nexis cannulated Ø2.7 drill bit for 4.0mm screws	1
<b>4</b>	XRE01007	Nexis <sup>®</sup> / PECA <sup>®</sup> -C Ø3.7 countersink for 4.0mm screws	1
<b>5</b>	XGA01004	Nexis graduated ruler 100mm	1