



AIRLOCK® FOOT TRAY OVERHEAD GUIDE



UNIVERSAL INSTRUMENTS

#	DESCRIPTION	PART NO.	QTY
-	TRAY	ACC1016P0001	1
-	LID	ACC1016P0003	1
1	FOOT SUPPORT PLATE	ACC1002P0008	1
2	AO HANDLE	XHA01001	1
3	DEPTH GAUGE	XGA01002	1
4	K-WIRE HOLDER	ACC1016P0004	1
-	K-WIRE Ø1,4 LG 150 TR/RD(1)	-	5
-	K-WIRE Ø1,6 LG 180 TR/RD(2)	-	5
5	CLEANING PIN Ø1,4	XKW01002	1

¹⁾K-wire supplied separately - Medetechnik® K-wire (33-T10-R-14-150) or Novastep® K-wire (CKW01010) are available depending on your market.

²⁾K-wire supplied separately - Medetechnik® K-wire (33-T10-R-16-180) or Novastep® K-wire (CKW01004) are available depending on your market.

AIRLOCK® INSTRUMENTATION

#	DESCRIPTION	PART NO.	QTY
6	T8 AO SCREWDRIVER TIP	XSD02002	2
7	SPHERIC POSITIONING PIN	XPP01003	2
8	SCREW LENGTH INDICATOR	XGA01003	1
9	PRESSLOCK® - DRILL GUIDE FOR SCREW Ø 3.5	XDG01023	1
10	OBLONG DRILL GUIDE FOR Ø3 SCREW	XDG01022	2
11	LOCKING DRILL GUIDE FOR SCREW Ø3 AND Ø3.5	XDG01021	1
12	DRILL BIT Ø2	XDB01021	2
13	DRILL BIT Ø2.5	XDB01022	OPTION
14	DRILL GUIDE HOLDER	XMS01004	1
15	POLYAXIAL DRILL GUIDE FOR SCREW Ø3 AND Ø3.5	XDG01020	1

REAMERS & PLATE BENDERS

#	DESCRIPTION	PART NO.	QTY
16	PLATE BENDERS	XMS01010	2
17	CONVEX REAMER Ø18	XRE01010	1
18	CONVEX REAMER Ø20	XRE01005	1
19	CONVEX REAMER Ø22	XRE01006	1
20	CONCAVE REAMER Ø18	XRE01011	1
21	CONCAVE REAMER Ø20	XRE01003	1
22	CONCAVE REAMER Ø22	XRE01004	1

NEXIS® INSTRUMENTATION

#	DESCRIPTION	PART NO.	QTY
23	SNAP-OFF SCREWDRIVER TIP	XSD03001	1
24	T10 AO SCREWDRIVER TIP	XSD04001	1
25	COUNTERSINK	XRE01007	1
26	CANNULATED DRILL BIT Ø2.7	XDB01007	1
27	DOUBLE DRILL GUIDE FOR SCREWS Ø4	XDG01009	1
28	RULER LG 150	XGA01009	1
29	T20 AO SCREWDRIVER TIP	XSD05001	OPTION
30	COUNTERSINK Ø4.9	XRE01008	OPTION
31	CANNULATED DRILL BIT Ø3.2	XBD01009	OPTION
32	DOUBLE DRILL GUIDE FOR SCREWS Ø4	XDG01015	OPTION
-	RULER LG 180/200	XGA01007	OPTION

PECA®-C INSTRUMENTATION

#	DESCRIPTION	PART NO.	QTY
33	EXACT2-T15 AO SCREWDRIVER TIP	XSD09002	OPTION
-	AO DRILL BIT Ø3.2	XDB01030	OPTION
-	K-WIRE Ø1,6 LG 150 TR/RD ⁽³⁾	-	OPTION

⁽³⁾K-wire supplied separately - Medetechnik® K-wire (33-T10-R-16-150) or Novastep® K-wire (CKW01003) are available depending on your market.

MTP TRIAL PLATES

#	DESCRIPTION	PART NO.	QTY
34	TRIAL PLATE HOLDER	ACC1006P0009	1
-	MTP MEDIUM TRIAL PLATE LEFT	XTI01001	1
-	MTP MEDIUM TRIAL PLATE RIGHT	XTI01002	1
-	MTP SHORT TRIAL PLATE LEFT	XTI01301	1
-	MTP SHORT TRIAL PLATE RIGHT	XTI01302	1
-	MTP LONG TRIAL PLATE LEFT	XTI01401	1
-	MTP LONG TRIAL PLATE RIGHT	XTI01402	1

FUSION TRIAL PLATES

#	DESCRIPTION	PART NO.	QTY
35	TRIAL PLATE HOLDER	ACC1006P0007	1
-	FUSION STRAIGHT SHORT TRIAL PLATE	XTI04117	1
-	FUSION STRAIGHT MEDIUM TRIAL PLATE	XTI04120	1
-	FUSION STRAIGHT LONG TRIAL PLATE	XTI04123	1
-	FUSION H SMALL TRIAL PLATE	XTI04217	1
-	FUSION H MEDIUM TRIAL PLATE	XTI04220	1
-	FUSION H LONG TRIAL PLATE	XTI04223	1

LISFRANC TRIAL PLATES

#	DESCRIPTION	PART NO.	QTY
36	TRIAL PLATE HOLDER	ACC1006P0012	1
-	LISFRANC H SMALL TRIAL PLATE	XTI05010	1
-	LISFRANC H MEDIUM TRIAL PLATE	XTI05020	1
-	LISFRANC H LARGE TRIAL PLATE	XTI05030	1
-	LISFRANC T SHORT TRIAL PLATE	XTI05040	1
-	LISFRANC T LONG TRIAL PLATE	XTI05050	1

UTILITY TRIAL PLATES


#	DESCRIPTION	PART NO.	QTY
38	TRIAL PLATE HOLDER	ACC1006P0010	1
-	UTILITY 2 HOLES TRIAL PLATE	XTI04016	1
-	UTILITY 3 HOLES TRIAL PLATE	XTI04022	1
-	UTILITY 4 HOLES TRIAL PLATE	XTI04028	1
-	UTILITY 5 HOLES TRIAL PLATE	XTI04034	1
-	UTILITY 6 HOLES TRIAL PLATE	XTI04040	1
-	UTILITY 7 HOLES TRIAL PLATE	XTI04046	1

LAPIDUS TRIAL PLATES

#	DESCRIPTION	PART NO.	QTY
37	TRIAL PLATE HOLDER	ACC1006P0011	1
-	LAPIDUS SHORT TRIAL PLATE RIGHT	XTI03301	1
-	LAPIDUS SHORT TRIAL PLATE LEFT	XTI03302	1
-	LAPIDUS LONG TRIAL PLATE RIGHT	XTI03401	1
-	LAPIDUS LONG TRIAL PLATE LEFT	XTI03402	1
-	PLANTAR LAPIDUS SHORT TRIAL PLATE RIGHT	XTI08101	1
-	PLANTAR LAPIDUS SHORT TRIAL PLATE LEFT	XTI08102	1
-	PLANTAR LAPIDUS LONG TRIAL PLATE RIGHT	XTI08201	1
-	PLANTAR LAPIDUS LONG TRIAL PLATE LEFT	XTI08202	1

MT BASE TRIAL PLATES - OPTION

#	DESCRIPTION	PART NO.	QTY
39	TRIAL PLATE HOLDER	ACC1006P0008	1
-	MT CLOSING WEDGE TRIAL PLATE LEFT	XTI02010	1
-	MT OPEN WEDGE 3 TRIAL PLATE LEFT	XTI02013	1
-	MT OPEN WEDGE 4 TRIAL PLATE LEFT	XTI02014	1
-	MT OPEN WEDGE 5 TRIAL PLATE LEFT	XTI02015	1
-	MT CLOSING WEDGE TRIAL PLATE RIGHT	XTI02020	1
-	MT OPEN WEDGE 3 TRIAL PLATE RIGHT	XTI02023	1
-	MT OPEN WEDGE 4 TRIAL PLATE RIGHT	XTI02024	1
-	MT OPEN WEDGE 5 TRIAL PLATE RIGHT	XTI02025	1



DISTRACTOR & COMPRESSOR TRAY - OPTIONAL

#	DESCRIPTION	PART NO.	QTY
-	DISTRACTOR TRAY	ACC1016P0002	1
40	DISTRACTOR K-WIRES HOLDER	ACC1016P0005	1
-	THREADED K-WIRE Ø2,5 - OPENING WEDGE DISTRACTOR - STERILE ⁽¹⁾	SKW04001	4
-	THREADED K-WIRE Ø1,6 LG 140 TR/RD - STERILE ⁽¹⁾	SKW05003	5
-	THREADED K-WIRE Ø2,5 LG 140 TR/RD - STERILE ⁽¹⁾	SKW05004	5
41	CLOSED ARMS DISTRACTOR	XFP01006	1
42	OPEN ARMS DISTRACTOR	XFP01008	1
43	OPENING WEDGE DISTRACTOR	XFP01011	1
44	COMPRESSIVE FORCEPS	XFP01012	1

⁽¹⁾Threaded k-wire supplied separately.

enovis

T +33 (0) 2 99 33 86 50 F + 33 (0) 9 70 29 18 95

Legal manufacturer: Novastep® S.A.S
 2 Allée Jacques Frimot | 35000 Rennes | France
 contact-intfa@enovis.com
 www.int.novastep.life
 Copyright © 2024 Enovis Foot and Ankle

Carefully read the enclosed Instructions For Use (IFU) and all packaging label information. Devices: Implants: Class IIb-CE1639 / Instruments: Class I / Class Ir-CE1639 / Class IIa-CE1639.

REFERENCE: AIR-TL02-ED4-02-26-EN