

# NEXIS® MIS TRAY OVERHEAD GUIDE



## NEXIS® MIS INSTRUMENTATION

| # | DESCRIPTION                             | PART NO.     | QTY |
|---|-----------------------------------------|--------------|-----|
| 1 | TRAY                                    | ACC1018P0001 | 1   |
| 2 | LID                                     | ACC1018P0002 | 1   |
| 3 | EXACT2-T8 SCREWDRIVER TIP               | XSD02006     | 2   |
| 4 | COUNTERSINK Ø2.8                        | XRE01024     | 1   |
| 5 | AO HANDLE                               | XHA01001     | 1   |
| 6 | FINE SURGICAL HANDLE <sup>(1)</sup>     | -            | 1   |
| 7 | PERIOESTAL ELEVATOR SINGLE TIP          | XMS01011     | 1   |
| 8 | RULER LG 100/150                        | XGA01013     | 1   |
| 9 | K-WIRE Ø1.2 LG 100 TR/RD <sup>(2)</sup> | -            | 5   |

<sup>(1)</sup> Reference SF13 supplied separately - availability depending on your market.

<sup>(2)</sup> K-wire supplied separately - Medetechnik® K-wire (33-T10-R-12-100) or Novastep® K-wire (CKW01014) are available depending on your market.

## OPTIONAL - NEXIS® MIS INSTRUMENTATION


| #  | DESCRIPTION                             | PART NO. | QTY |
|----|-----------------------------------------|----------|-----|
| 10 | REMOVAL EXACT2-T8 AO SCREWDRIVER TIP    | XSD02007 | 1   |
| 11 | CANNULATED DRILL BIT Ø1.9               | XDB01025 | 1   |
| 12 | K-WIRE Ø1.2 LG 150 TR/RD <sup>(3)</sup> | -        | 5   |

<sup>(3)</sup> K-wire supplied separately - Medetechnik® K-wire (33-T10-R-12-150) or Novastep® K-wire (CKW01015) are available depending on your market.

## OPTIONAL - NEXIS® Ø2,3 INSTRUMENTATION

| #  | DESCRIPTION                             | PART NO. | QTY |
|----|-----------------------------------------|----------|-----|
| 13 | T7 AO SCREWDRIVER TIP                   | XSD01001 | 2   |
| 14 | K-WIRE Ø1.2 LG 150 TR/RD <sup>(4)</sup> | -        | 5   |

<sup>(4)</sup> K-wire supplied separately - Medetechnik® K-wire (33-T10-R-09-100) or Novastep® K-wire (CKW01008) are available depending on your market.



T +33 (0) 2 99 33 86 50 F +33 (0) 9 70 29 18 95

Legal manufacturer: Novastep® S.A.S  
2 Allée Jacques Frimot | 35000 Rennes | France  
contact-intfa@enovis.com  
www.int.novastep.life  
Copyright © 2024 Enovis Foot and Ankle

Carefully read the enclosed Instructions For Use (IFU) and all packaging label information. Devices: Implants: Class IIb-CE1639 / Instruments: Class I / Class Ir-CE1639 / Class IIa-CE1639.

REFERENCE: NEX-MIS-TL-ED4-06-25-EN