



Brochure

-

Sterile Percutaneous Burrs



Keys & Features

The Novastep Percutaneous Burrs are designed specifically for Mini Invasive Surgery needs. Cylindrical and conical burrs are made of hardened surgical grade stainless steel, highly resistant to wear, providing an optimal cutting quality.



Hammertoe, Akinette:

Shannon Corta Ø 2.0 x 8 mm



Akin, DMMO:

Shannon Recta Ø 2.0 x 12 mm
Shannon Helical Ø 2.0 x 12 mm



Bunion, Joint prep:

Shannon Longa Ø 2.2 x 22 mm



Calcaneal Slide:

Shannon Larga Ø 3.0 x 20 mm (7cm)
Shannon X-Larga Ø 3.0 x 30 mm (10cm)



Cheilectomy, Osteophytes:


Wedge Ø 3.1 x 13 mm



Cheilectomy, Osteophytes:

Wedge Ø 4.1 x 13 mm

References



Designation	Reference
Shannon Corta Ø2.0 x 8	CRE12008
Shannon Recta Ø2.0 x 12	CRE12012
Shannon helical Ø2.0 x 12	CRE12212
Shannon Longa Ø2.2 x 22	CRE12222
Shannon Larga Ø3.0 x 20	CRE13020
Shannon X-Larga Ø3.0 x 30	CRE13030
Wedge Ø3.1	CRE23113
Wedge Ø4.1	CRE24113

Recommandation :

Carefully read the enclosed Instructions For Use (IFU) and all packaging label information. Devices: Instruments: Class IIa-CE1639.

Novastep:

2, Allée Jacques Frimot
35000 RENNES France
Tel : + 33 (0) 2 99 33 86 50 / Fax : + 33 (0) 9 70 29 18 95
contact@novastep-ortho.com / www.novastep-ortho.com

Per-L-Ed1-11-21-EN