

percutaneous burrs



BROCHURE

FOREFOOT
MIDFOOT
REARFOOT

Sterile percutaneous burrs



*Creating
Better
Together™*

percutaneous burrs

Keys & features

The Novastep Percutaneous Burrs are designed specifically for Mini Invasive Surgery needs. Cylindrical and conical burrs are made of hardened surgical grade stainless steel, highly resistant to wear, providing an optimal cutting quality.



Hammertoe, Akinette:
Shannon Corta Ø 2 Lg 8 mm



Akin, DMMO:
Shannon Recta Ø 2 Lg 12 mm
Shannon Helical Ø 2 Lg 12 mm



Bunion, Joint prep:
Shannon Longa Ø 2.2 Lg 22 mm



Calcaneal Slide:
Shannon Larga Ø 3 Lg 20 mm (7cm)
Shannon X-Larga Ø 3 Lg 30 mm (10cm)




Cheilectomy, Osteophytes:
Wedge Ø 3.1 Lg 13 mm



Cheilectomy, Osteophytes:
Wedge Ø 4.1 Lg 13 mm

References

Ref	Description
CRE12008	Shannon Corta Ø 2 Lg 8 mm
CRE12012	Shannon Recta Ø 2 Lg 12 mm
CRE12212	Shannon Helical Ø 2 Lg 12 mm
CRE12222	Shannon Longa Ø 2.2 Lg 22 mm
CRE13020	Shannon Larga Ø 3 Lg 20 mm
CRE13030	Shannon X-Larga Ø 3 Lg 30 mm
CRE23113	Wedge Ø 3.1 Lg 13 mm
CRE24113	Wedge Ø 4.1 Lg 13 mm



Please note:

Carefully read the enclosed Instructions For Use (IFU) and all packaging label information. Devices: Implants: Class IIb-CE1639 / Instruments: Class I / Class Ir-CE1639 / Class IIa-CE1639.

Novastep:

2, Allée Jacques Frimot - 35000 RENNES - France
Tel. : + 33 (0) 2 99 33 86 50 / Fax : + 33 (0) 9 70 29 18 95
contact@novastep-ortho.com / www.int.novastep.life

Reference: Per-L-Ed2-07-24-EN