

EXray[®] 20

Snap-off
Threaded pin



SynchroMedical

The interphalangeal reference

distributed by **enovis** FOOT & ANKLE

EXray[®] 20 PEEK

Snap-off threaded pin

DESCRIPTION

EXRAY[®] 20 snap-off threaded pin has been developed as a resectable peek pin. The seat of the implant is prepared with a 2 mm drill. The pin consists of a trocar tip, a threaded portion and a portion with a narrower diameter that can be connected to a power tool.

The EXRAY[®] 20 has a non-compressive \varnothing 2.5mm continuous thread. The 3 markings are positioned at 15 mm, 20 mm and 25 mm. The full length of the thread is 30 mm. Once implanted, the excess length is removed with a saw blade.

EXRAY[®] 20 snap-off threaded pins are manufactured in PEEK as per standard ASTM F2026.

- Simple preparation
- Variable screw length
- Used to support stabilization of osteotomies
- Simplified revision
- PEEK threaded pins are radiolucent
- The elasticity modulus is close to the one of bone



ORDERING INFORMATION

Implant

PW2060.....L 60 mm



Instruments

	AO20ST100	DRILL
	PW9010	PROTECTIVE BARREL
	HV9016	SCREW EXTRACTOR

The counterclock screw extractor (Ref. HV9016) is used for revision

All rights reserved © Pictures non-contractual